



Tackle Your Vision

Lucy Bronze, MBE, chose the EVO ICL procedure to help her get rid of the daily hassles of contact lenses and glasses, for the long term.



The visual freedom you deserve.

Ditch foggy glasses and dry contacts.

The EVO ICL lens allows you to unlock visual freedom and a life rich with new experiences.

EVO can help you **see the world the way it's meant to be seen** – with sharp, clear vision.¹

A proven history
of ICL success...



Over

3,000,000+

ICLs distributed
worldwide.



Exceptional Patient Satisfaction

99.4%

of patients surveyed say they
would have the procedure again.²

ICL has over

30 years

of experience
in the eye


Benefits of EVO ICL™

Exceptional Vision Day & Night



Experience what life has to offer throughout the day **without compromising your night vision.**^{1,3}

Day or night—you can rely on EVO for exceptional vision.^{1,3}



.....| **In Harmony with Your Natural Eye**



Similar to adding a lens to a camera, EVO is additive and works together with the natural parts of your eye, to improve your vision.

.....| **No Dry Eye Syndrome**



Clinically proven to not cause or worsen dry eye syndrome.⁴

.....| **Your Future's Bright & Protected**



EVO is made of Collamer®, a material shown to be safe in the eye. Collamer provides UV protection and contains collagen, a naturally occurring protein in the body.

Created for your well-being

Made with Collamer — a safe, biocompatible and premium material

Collamer is biocompatible with your eye, and is made with collagen, **a protein which naturally occurs in your body**. Shown to be safe and effective, and produced with high quality materials, Collamer is used exclusively by STAAR Surgical in EVO.

..... | **Soft & Flexible**



EVO is ideal for your eye because of its **soft and flexible** structure. You won't even know it's there, apart from the fact that you have great vision.

Keep your peace of mind

.....| A Removable Option That Gives You Flexibility For The Future




While EVO is designed to permanently correct your vision, **the lens can be removed** by a doctor, if necessary, or to keep pace with advancing technology and your future vision needs.



How EVO ICL™ works


EVO functions in a similar way to glasses and contact lenses to correct myopia and myopia with astigmatism. **EVO is placed behind the iris (the colored part of the eye) and right in front of the natural lens, helping you to achieve sharp, clear vision.¹**

Nearsightedness (Myopia)




Light rays focus in front of the retina

Astigmatism




Light rays fail to focus at a single point on the retina

After Treatment



EVO ICL

EVO focuses light onto the retina creating clear vision



Quick and Safe Procedure

The EVO procedure is as quick or faster than watching an episode of your favorite binge-worthy series and has a short recovery time.

Getting ready for visual freedom!


Before the day of your EVO procedure, your doctor will perform a series of standard tests to measure your eye's unique characteristics.

Time for your EVO lens!

When you arrive for your procedure, your doctor will administer eye drops to dilate and numb your eye.


The Procedure Follows Three Easy Steps:

1




Your doctor will create a small opening at the edge of your cornea to insert EVO.

2



EVO is folded and inserted through this small opening.

3



Once the lens is inserted, your doctor will make any necessary adjustments to ensure proper positioning in the eye.

EVO is not visible to you or others.

The procedure is complete!

Many patients notice improved vision nearly immediately.

After-care: Consult with your doctor for further information on the recovery process.

Why Lucy Bronze chose EVO

Lucy Bronze' struggle with poor vision

Lucy Bronze chose the EVO ICL procedure to correct her vision, which has been game changing for her, both on and off the pitch. As a professional footballer, she must react immediately during matches and having sharp, clear vision is essential.

Previously, she struggled with short-sightedness, wearing glasses or contact lenses from the age of 9. Constantly breaking her glasses or experiencing issues with contact lenses, such as irritation and dry eyes, caused her much frustration.

With EVO ICL, she no longer has to worry about these problems. Her biggest regret is not undergoing the procedure sooner, as it has given her the freedom to play and live her life confidently, without vision-related concerns.





“

The EVO ICL procedure was even better than expected, it was pain free and super quick for me. My recovery was unbelievable, I was back playing sport, training and playing football within 2-3 weeks.



Lucy Bronze, MBE
Professional Footballer



IMAGINE YOUR LIFE WITH VISUAL FREEDOM

DISCOVER MORE ABOUT EVO ICL TODAY!



Important Safety Information For EVO/EVO+ ICL

The EVO/EVO+ ICL is designed for the correction/reduction of up to -20 diopters (D) of nearsightedness with up to 6 D of astigmatism for patients who are 21 to 60 years of age. Implantation of the EVO/EVO+ ICL is a surgical procedure, and as such, carries potentially serious risks. The following represent potential complications/adverse events: additional surgeries, cataract formation, transient or persistent loss of best corrected vision, raised pressure inside the eye, loss of cells on the innermost surface of the cornea, conjunctival irritation, corneal swelling, conjunctival irritation, endophthalmitis (total eye infection), significant glare and/or halos around lights, hyphema (blood in the eye), hypopyon (pus in the eye), eye infection, EVO/EVO+ ICL dislocation, macular edema, non-reactive pupil, pupillary block glaucoma, severe inflammation of the eye, iritis, uveitis, vitreous loss and corneal transplant. Before considering EVO/EVO+ ICL surgery you should have a complete eye examination and talk with your eye care professional about the EVO/EVO+ ICL procedure, especially the potential benefits, risks, and complications. You should discuss the time needed for healing after surgery. For additional information with potential benefits, risks and complications please visit DiscoverICL.com

1. Martínez-Plaza E, López-Miguel A, López-De La Rosa A, et al. Effect of the EVO+ Visian Phakic Implantable Collamer Lens on Visual Performance and Quality of Vision and Life, *Am J Ophthalmol* 2021;226: 117–125. **2.** Packer M. The Implantable Collamer Lens with a central port: review of the literature. *Clin Ophthalmol.* 2018;12:2427-38. **3.** Parkhurst GD. A prospective comparison of phakic collamer lenses and wavefront-optimized laser-assisted in situ keratomileusis for correction of myopia. *Clin Ophthalmol.* 2016 Jun 29;10:1209-15. **4.** Naves, J. Carracedo, G. Cacho-Babillo, I. Diadenosine Nucleotid Measurements as Dry-Eye Score in Patients After LASIK and ICL Surgery. Presented at American Society of Cataract and Refractive Surgery (ASCRS) 2012.

staar.com

© 2024 STAAR® Surgical Company

For more information, visit: **evoicl.com**

UK-EVO ICL-24-0059

STAAR® Surgical Company
25651 Atlantic Ocean Drive
Lake Forest, CA 92630
(800) 352-7842